

ThevoSiiS – 2018 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 | Fax 866.870.0801

Date	Buyer	PO#	Mark For
Address			
Contact	Phone ()	Fax ()	

Step 1 – Seat Sizes



Seat: - firm seat block provides maximum sensory stimulation, - elongated shape of the seat helps the forward facilitation of the arms. By bracing the arms on the set block an interim impact change for the spine can be achieved, - seat block forces the user into a physiological seating position, - seat block is sloped in the back part in order to slightly tilt the pelvis.

Back: - small and firm iliac crest support with slight lean against the body helps the stabilisation, - iliac spine support and impulse generator are adjustable, - point-shaped impulse transmitters give an orientation for the back position and stimulate the erection and orientation process.

Base frame: - adjustable in height, - lever brake and height adjustment can be operated independently by the user, - therapy chair is tilt resistant, - push bar enables an assistant to also carry the chair.

Areas of application: cerebral apraxia, cognitive processing deficiencies/ ADD, tonus dysregulation deficiencies.

HCPCS	£ mini	£ size 2	£ narrow	£ plus	£ high
ThevoSiiS E1399	8301	8305	8308	8306	8307
	\$1749	\$1898	\$1898	\$1898	\$2049



ThevoSiiS	mini	size 2	narrow	plus	high
Item code	8301	8305	8308	8306	8307
Seat height	12.6 - 15"	16.5 - 24"	16.5 - 24"	19.3 - 28"	20.9 - 31.5"
Seat width	5.9"	7.5"	6.3"	8.7"	7.5"
Impulse generators					
Height adjustment	10.8"	7.1"	7.1"	7.1"	10.8"
Depth adjustment	1.4"	-	-	-	1.4"
Iliac crest support					
Height adjustment	3.2"	3.2"	3.2"	3.2"	3.2"
Depth adjustment	1.4"	1.4"	1.4"	1.4"	1.4"
Total dimensions (wxhxl)	21.3 x 19.7 x 26.8"	23.6 x 34.7 x 22"	23.6 x 34.7 x 22"	23.6 x 34.7 x 22"	23.6 x 34.7 x 22"
Push bar height		34.7"	34.7"	34.7"	34.7"
Wheel size	3.9"	3.9"	3.9"	3.9"	3.9"
Wheel distance front	20.9"	22.8"	22.8"	22.8"	22.8"
Wheel distance rear	21.3"	21.7"	21.7"	21.7"	21.7"
Max. weight capacity	49 lbs	265 lbs	265 lbs	265 lbs	265 lbs
Weight	24 lbs	31 lbs	31 lbs	31 lbs	33 lbs
Body size	max. 49.2"	max. 53.2"	max. 53.2"	max. 55.1"	max. 57.1"

Step 2 – General Options

General Accessories	HCPCS	Unit	Price
Front wheels, lockable*	K0108	Each	£ 8321 \$95
Impulse generators, depth adjustable*	K0108	Each	£ 8320 \$30

*items sold individually, for new orders two of this item is required.