

**ThevoTherapy – 2021 US Retail prices** info@thomashilfen.com | Phone 866.870.2122 | Fax 866.870.0801

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

**Step 1 – Seat Sizes (Required options are indicated below by † symbol)**

	£ ThevoTherapy seat size 1		NA ThevoTherapy seat size 2		HCPCS	
†Seat	840011	\$365	845011	\$415	E2292	Individual body adapted, integrated seat/back with MIS-Microstimulation®; 122° recline gas spring operated; growth adjustable seat upholstery; articulated footrest height and angle adjustable; height, depth angle angle adjustable swing away head rest bracket; front facing seat; ventilated/washable upholstery; Dynamic Seating in 3rd Dimension.
†Back	840021	\$475	845021	\$525	E2293	
†Recline	840031	\$350	845031	\$350	E1014	
†Hip sup.	840032	\$200	845032	\$200	E0956	
†Calf panel	6635	\$35	6635	\$35	E0995	
†Footplate, flip up	840035	\$70	845035	\$70	K0040	
†Legrest, elevating	840036	\$200	845036	\$200	E0990	



Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Seat to Head	Max Height	Hip Angle
Size 1	6 ¾ – 12”	7 ¼ – 12”	7 ¼ – 12 ¼”	18 ½ – 25 ¾”	54”	90° - 122°
Size 2	10 – 14 ¼”	10 ¾ – 15”	9 ¾ – 14 ¼”	25 ¼ – 30 ¾”	63”	90° - 122°



**Step 2 – Frame Options**

**Hi-low chassis** - therapy frame, height/tilt easily adjustable by gas spring, 4 swiveling wheels, back wheels lockable

	HCPCS	Size 1		Size 2	
Hi-low indoor chassis, pedal controlled gas spring	E1233	£ 840041	\$1599	£ 845041	\$1599
Push bar	K0108	£ 8439	\$225	£ 8439	\$225

**Step 3 – Positioning & General Options**

Lower Body Positioning		HCPCS	Unit	Small	>>>Medium	>>>Large	Price
Foot straps with cam buckle		E0952	Pair	£ 6648			\$30
Ankle Huggers, Bodypoint®		K0108	Pair	£ 213151			\$205

**Mid Body Positioning**

Abduction block, padded		E0957	Each	£ 6620	£ 6720	\$138
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420	£ 8470	\$80
	4 point, padded, SR buckle	E0978	Each	£ 8421	£ 8471	\$130
Crotch strap (fits 8420,8470,8421,8471)		E0978	Each	£ 6341		\$28
Arm rest, pair, height, depth & angle adjustable		E0973	Pair	£ 8431	£ 8478	\$375
Tray, fold-away, black, height adjustable,		E0950	Each	£ 8425	£ 8475	\$315
Hip support, high, padded, l&r, adjustable		E0956	Pair	£ 8423	£ 8473	\$149

**Upper Body Positioning**

Chest harness supports	Butterfly style, neoprene **	E0960	Each	£ 6385	£ 6585	\$165
	H-style, posture support, AEL **	E0960	Each		£ 14813	\$185
	PivotFit shoulder harness, dynamic, Bodypoint® **	E0960	Each		£ 214270	\$185
	PivotFit shoulder harness, Bodypoint® **	E0960	Each		£ 214250	\$165
Trunk supports	Lateral, dynamic, padded, l&r, height adj.	E0956	Pair	£ 8422	£ 8472	\$255
	Rigid lateral padded, l&r, height adjustable	E0956	Pair	£ 8426	£ 8476	\$230
	Flexible lateral padded l&r, w/chest strap, adj.	E0956	Pair		£ 8427	\$285
Head support	Shell*	E0955	Each		£ 84240	\$255
	Swing-away supp. bracket shell*	E1028	Each		£ 84241	\$150
	Neck support*	E0955	Each		£ 84290	\$255
	Swing-away supp. bracket neck*	E1028	Each		£ 84291	\$150

**General Accessories**

Pelvis support, padded		K0108	Each	£ 8432	£ 8479	\$89
Belt bow (mounting bracket for harness supports) **		K0108	Each	£ 6371	£ 6571	\$119

\* Head rest with bracket is required option because seat back grows through head rest height adjustment  
 \*\* Belt bow is required for usage of chest harness supports

Technical Specs #		
ThevoTherapy	Size 1	Size 2
Capacity	77 lbs	88 lbs
Unit weight	44 lbs	47 lbs
Max seat tilt	31°	31°
Seat direction	front facing	
Upholstery	Trevira-CS	
Floor to handle	31-41”	
Floor to seat	10-26”	
Max width	21”	
Max length	29”	
Fold width	21”	
Fold height	28”	
Fold length	29”	
Brakes	foot operated	
Front wheels	3.9”	
Rear wheels	3.9”, lockable	

# More details in Technical Manual

