



Rehab-equipment
INCLUSIVE 

ACCESSORIES
(more available)



The practical padded swivel base, the Seatfix connection for ISOFIX and the footrest give additional stability.



The seat wedge inside and below provide a relaxed resting position during the journey. The abduction block offers extra positioning for children whose legs cross or "scissor".



The table is a support surface and an additional impact shield that also prevents the child from attempting to unbuckle the 5-point positioning harness.

Measurements & weights	RECARO Monza Nova 2 Reha
Seat depth	30 cm / 11.8" (43 cm / 16.9" with seat depth extension)
Seat width	front: 29 cm / 11.4" rear: 24 cm / 9.4" (20 cm / 7.9" with seat wedge, inside)
Back height	60 - 77 cm / 23.6 - 30.3"
Lower leg length <small>(footrest mounted on swivel plate)</small>	60 - 77 cm / 23.6 - 30.3"
Shoulder height (child)	33 - 48 cm / 13 - 18.9"
Overall dimensions (w x h x d)	54 x 66 x 42 cm / 21.3 x 26 x 16.5"
Weight of seat approx.	7 kg / 15.4 lb
Carrying capacity	15 - 50 kg / 33.1 - 110.2 lbs
Children height	94 - 150 cm / 37 - 59"



Products & accessories	item code
RECARO Monza Nova 2 Seatfix Reha	841
Abduction block	804
Swivel base with footrest adapter* (incl. adapter plate)	820/1
Recommended options for stability:	
• Seatfix connection (for swivel base) option 1	818
• Stabilising bow, rear (for swivel base) option 2	819
Stand, front (for swivel base)	817
Footrest, long* (depth: 18 cm / 7.1") (tiltable and adjustable in angle)	805
Footrest, short* (depth: 12 cm / 4.7") (tiltable and adjustable in angle)	806
Seat depth extension* (+ 13 cm / + 5.1")	808
Seat wedge, inside	809
Seat wedge, below* (+ 10°)	810
Table	811
Soft harness strap cover, grey (pair)	815

* Footrest and seat depth extension cannot be used in combination.



EXTRA SUPPORTIVE
for children with special needs

“ The safety of my child is most important to me. Generally, I am a careful driver, but you never know. In case of an accident I have to be 100% sure that the car seat protects my child against injuries.

”

Susanne Bessel, mother




RECARO Monza Nova 2 Reha is reliably tested many times

- ✓ FMVSS 213 (NHTSA Standard for USA Car Seats)
- ✓ UN/ECE R 44/04 (United Nations Economic Commission for Europe – regulation for child restraint systems)
- ✓ DEKRA (tests for „Special Needs Restraint“) incl. testing of all rehab accessories
- ✓ Manufacturer’s guarantee



Pay attention to this seal

Thus you can be sure that the **rehab seat** successfully passed the necessary overload-test up to 50 kg / 110.2 lbs according to ECE R 44/04.



RECARO Monza Nova 2 Reha
Lightweight booster seat exceeds highest safety standard



Quick and easy seat transfer from one vehicle to another.



Lightest rehab seat in its class that you can carry with one finger ... Yet it carries a child up to 15 to 50 kg / 33 to 110 lbs.



Improved strength and static-free side pads reduce crash-generated forces to the head, neck and shoulders by up to 30%.

30%

Extra-firm armrests for even greater protection in the event of side impact.