

# New rehabilitation stroller offers safety in comfort

## JAZZ EASyS combines safety, ergonomics and a feel-good factor

Child patients are on the move for several hours a day in their rehabilitation strollers. Whether accompanying their parents on shopping trips or on the way to play with their little friends, the stroller is their constant companion and trusty steed. A good reason, then, to continually improve existing models to provide the best possible support for handicapped children and their families in their day to day life. As Gunnar Thomas, owner and managing director of Thomashilfen explains "It is our philosophy to offer the utmost safety in our rehabilitation strollers", the number one priority for his company. This is why the company's latest model of stroller, a practical combination of seat shell and chassis, also features the Balance-Safety-System: when adjusting the backrest, the whole seating unit is pushed forwards, thus stabilising the stroller's centre of gravity, so the BSS® offers considerable safety reserves, especially for larger children

### Smart combinations

With this new rehabilitation stroller, called the JAZZ EASyS, the motto "safety first" alone was far from sufficient for the manufacturer – here, it was also about smart combinations: many years of experience in the area of child rehabilitation and clever new ideas have been put into practice using innovative technology and the most up-to-date materials. The results speak

for themselves. "We have succeeded in combining optimum safety and maximum comfort" - Gunnar Thomas is proud of this model, which conforms to all international safety standards, from crash tests to fire-resistant properties. Thomas describes the advantages of this modern-designed stroller model "Children feel safe and sound in the stable and comfortable stroller, and parents soon notice how easy and practical it is to use", and lists the most important new features: "The adjustment angles in the hip and knee areas have been physiologically optimised and now run parallel to the movement axis of the hip and knee joint. This guarantees that supporting accessories can always be correctly installed for therapy purposes. This new feature is unique in the global market." As well as the ability to adjust the seat angle to take the strain off the spine and circulation, the seating system, which grows with your child, also allows you to flexibly adjust the angle of the backrest for an optimum sitting posture at all times.

### Feeling good made easy

Gunnar Thomas explains a further advantage of the new JAZZ EASyS: "To protect the kids from excessive sweating or chafing, even after long periods in the stroller, we have used an especially breathable spacer fabric in the cover of the soft padded seat unit". Another new feature is the seat's lockable suspension. The suspension, which can be locked when required, gently cushions against bumps and uneven surfaces. However, when making improvements in comfort, we considered not just the kids, but also the parents: "EASyS Klick-System" is the name of the system that allows you to securely attach the seat unit to the adapter or to quickly release it, at the



*The JAZZ EASyS stroller has also an appealing classical design.*

press of a button. Once "unclicked", the stroller, which weighs only around 15 kilos and yet is extremely robust, will fit in any car boot.

### Testing and customising

As with all rehabilitation strollers in the EASyS range, this new model can also be quickly and easily adapted to the needs of your child. Our company offers help and advice in finding the ideal solution for your individual needs. As a special service, the manufacturer offers the opportunity to test this modern rehabilitation stroller in day to day life, free of charge and without obligation, letting parents and children see for themselves how enjoyable everyday life can be with the JAZZ EASyS.

*Further information on the JAZZ EASyS rehabilitation buggy is available directly from the manufacturer Thomashilfen, Walkmühlenstr. 1, D - 27432 Bremervörde or at the number ++49 4761 / 8860*



*Every day can be a good experience, when children feel confident in their stroller.*