

EASyS Advantage – 2022 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 | Fax 866.870.0801

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat Sizes (Required options are indicated below by † symbol)

EASyS seat EASyS seat



	size 1	size 2	HCPCS	Description
†Seat	660961 \$420	670961 \$420	E2292	Solid contoured, integrated seat/back, 180° recline, growth adjustable; seat depth/ width adjustable from underneath, articulated foot plate height and angle adjustable; knee angfront/rear facing seat, easy to change; ventilated/washable upholstery ; knee angle adjustment with release lever; cozy comfortable seat cushion, BSS-Balance Safety System; calf strap; 5 different color options
†Back	660962 \$518	670962 \$518	E2293	
†Recline	660963 \$415	670963 \$415	E1014	
†Hip sup.	660964 \$230	670964 \$230	E0956	
†Calf panel	6635 \$45	6735 \$45	E0995	
†Foot plate, flip up	660965 \$90	670965 \$90	K0040	
†Legrest, elevating	660966 \$200	670966 \$200	E0990	

Seat Size Chart	Lower Leg	Seat Depth	Seat width	Seat to Head	Max Height	Hip Angle
Size 1	7 – 11½"	7½ – 11¾"	7½ – 11¾"	16½ – 26¾"	49"	90° - 180°
Size 2	8½ – 14½"	10½ – 15"	10 – 14½"	21½ – 29"	60"	90° - 180°

†Upholstery colors	Red/ Grey	Blue/ Grey	Green/ Grey	Magenta/ Grey	Grey/ Grey
Size 1	<input type="checkbox"/> 6247	<input type="checkbox"/> 6248	<input type="checkbox"/> 6249	<input type="checkbox"/> 6246	<input type="checkbox"/> 6245
Size 2	<input type="checkbox"/> 6447	<input type="checkbox"/> 6448	<input type="checkbox"/> 6449	<input type="checkbox"/> 6446	<input type="checkbox"/> 6445

Step 2a – Outdoor Frame Option (Required options are indicated below by † symbol)

Advantage-chassis - Folding aluminum frame, tilt-in-space, back frame suspension, face-to-face option, 0° safety lock out, cargo basket, unfold lock



	HCPCS	SIZE 1	SIZE 2
Advantage-chassis	E1234	<input type="checkbox"/> 62023 \$2295	<input type="checkbox"/> 64023 \$2295
†Castor pin lock, pair	K0073	67830 \$96	67830 \$96
†Bus tie-down loop kit (*)	K0108	67839 \$205	67839 \$205
Pneumatic wheels, full set	K0070	<input type="checkbox"/> 66461 \$249	<input type="checkbox"/> 66461 \$249
Frame padding, black	K0108	<input type="checkbox"/> 6622 \$115	<input type="checkbox"/> 6622 \$115
Frame padding, red/ black	K0108	<input type="checkbox"/> 66221 \$115	<input type="checkbox"/> 66221 \$115
Tip assist/ curb tipper	K0108	<input type="checkbox"/> 67502 \$89	<input type="checkbox"/> 67502 \$89
Under seat storage basket	K0108	NC	NC

EASyS Advantage Technical Specs	SIZE 1	SIZE 2
Seat capacity	77.2 lbs	88.2 lbs
Frame capacity	110.2 lbs	121.3 lbs
Unit weight	41 lbs	48 lbs
Max seat tilt	-10° to +35°	
Frame suspension	Yes	
Upholstery	Trevira-CS	
Foot angle	78° - 100°	
Size of foot rest	11" x 7.1"	12.6" x 8.3"
Floor to handle	28.2 – 46.1"	28.2 – 46.1"
Floor to seat	19"	19.5"
Max width	23.4"	27.2"
Fold length	33.5"	34"
Fold height	23.6"	27.2"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7.1" solid	
Rear wheels	10" solid	



* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Step 2b – Indoor Frame Option (Required options are indicated below by † symbol)

Q-chassis - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas springs, lockable rear wheels



	HCPCS	SIZE 1 & 2	
Q-chassis indoor hi/low, pedal controlled	E1399	<input type="checkbox"/> 6666	\$1999
†EASyS Adapter for Q-chassis	K0108	<input type="checkbox"/> 62951	\$130
Push bar for Q-chassis	K0108	<input type="checkbox"/> 6660	\$298

Q-chassis Technical Specs **

Capacity	88 lbs
Seat tilt	-7 to +31°
Floor to seat	11 - 28"
Seat direction	front facing

Step 3 – Positioning & General Options (Required options are indicated below by † symbol)

Lower Body Positioning	HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle	E0952	Pair		<input type="checkbox"/> 6648		\$35
Ankle Huggers, Bodypoint®	K0108	Pair		<input type="checkbox"/> 213151		\$205
Footrest pad, height 1¾"	K0108	Each	<input type="checkbox"/> 10622		<input type="checkbox"/> 10023	\$99

Mid Body Positioning						
Abduction block, padded	E0957	Each	<input type="checkbox"/> 6620		<input type="checkbox"/> 6720	\$160
Hip supports	2 point, padded, SR buckle	E0978	Each	<input type="checkbox"/> 8420	<input type="checkbox"/> 8470	\$95
	4 point, padded, SR buckle	E0978	Each	<input type="checkbox"/> 8421	<input type="checkbox"/> 8471	\$135
	Soft pelvic harness/ leg harness	E0978	Each	<input type="checkbox"/> 6614	<input type="checkbox"/> 6714	\$175
Crotch strap (fits 6629,8420,8470,8421,8471)	E0978	Each		<input type="checkbox"/> 6341		\$33
Tray, height adjustable, clear with picture window	E0950	Each	<input type="checkbox"/> 66281		<input type="checkbox"/> 67281	\$190
Mounting hardware tray - †required when ordering 6628/6728	E1028	Each	<input type="checkbox"/> 66282		<input type="checkbox"/> 67282	\$225
Arm rest, height, depth and angle adjustable	E0973	Pair	<input type="checkbox"/> 8431		<input type="checkbox"/> 8478	\$375

Upper Body Positioning							
Chest harness supports	Butterfly style, neoprene	E0960	Each	<input type="checkbox"/> 6385		<input type="checkbox"/> 6585	\$189
	Vest style	E0960	Each	<input type="checkbox"/> 6631	<input type="checkbox"/> 6632	<input type="checkbox"/> 6731	\$189
	H-style, posture support, small, Bodypoint®	E0960	Each		<input type="checkbox"/> 214806		\$199
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		<input type="checkbox"/> 214270		\$199
	PivotFit shoulder harness, Bodypoint®	E0960	Each		<input type="checkbox"/> 214250		\$189
Trunk supports	Rigid lateral, left & right, height adjustable	E0956	Pair	<input type="checkbox"/> 6634		<input type="checkbox"/> 6734	\$245
	Flexible lateral, l&r w/ chest strap, height adj.	E0956	Pair		<input type="checkbox"/> 6633		\$295
Head supports	Contoured pillow	E0955	Each		<input type="checkbox"/> 6824		\$152
	Rigid lateral, left & right, height adjustable	E0955	Pair	<input type="checkbox"/> 6624		<input type="checkbox"/> 6724	\$215

General Accessories							
Grab bar, padded, black	K0108	Each	<input type="checkbox"/> 6621		<input type="checkbox"/> 6721	\$125	
Footrest tube upholstery, grey	K0108	Pair		<input type="checkbox"/> 6230		\$28	
Vent tray	E1029	Each	<input type="checkbox"/> 6294		<input type="checkbox"/> 6494	\$450	
Medical necessities pack, black	K0108	Each		<input type="checkbox"/> 6639		\$80	
Oxygen tank holder, E-tank	E2208	Each		<input type="checkbox"/> 665012		\$210	
Straps	Secure folded seat/frame	K0108	Each		<input type="checkbox"/> 667011	\$30	
	Wrist "runaway" tether	K0108	Each		<input type="checkbox"/> 667111	\$30	
Weather protection	Canopy, sun/rain with window, black	K0108	Each		<input type="checkbox"/> 663761	\$325	
	Lower body windbreaker, black for 663761	K0108	Each		<input type="checkbox"/> 661852	\$99	
	Canopy, rain only, clear	K0108	Each		<input type="checkbox"/> 6638	\$245	
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		<input type="checkbox"/> 6739	\$125	
	Rain cape, grey	K0108	Each	<input type="checkbox"/> 6340		<input type="checkbox"/> 6540	\$163
	Sleeping bag, summer weight, black	K0108	Each	<input type="checkbox"/> 6619		<input type="checkbox"/> 6719	\$299
	Sleeping bag, winter weight, black	K0108	Each	<input type="checkbox"/> 6618		<input type="checkbox"/> 6718	\$299
Umbrella, angle adjustable	K0108	Each		<input type="checkbox"/> 6594		\$95	

*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement. ExoMotion LLC dba Thomashilfen | www.thomashilfen.us Effective Date January 1st 2022