

# EASyS Modular S – 2018 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 | Fax 866.870.0801

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

## Step 1 – Seat Sizes (Required options are indicated below by † symbol)

£ Modular S seat   £ Modular S seat



	size 1	size 2	HCPCS	Description
†Seat	661011	671011	\$475 \$525	E2294
†Back	661012	671012	\$475 \$525	E2293
†Recline	661013	671013	\$350 \$350	E1014
†Hip sup.	661014	671014	\$200 \$200	E0956
†Bracket	661015	671015	\$125 \$125	E1028
†Calf strap	6635	6735	\$35 \$35	E0995
†Foot plate	660965	670965	\$60 \$60	K0040

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Back Height#	Max Height	Hip Angle
Size 1	7.3 – 12"	7.5 – 11.8"	7.5 – 10.2"	11 – 16.5"	50"	90° - 140°
Size 2	7.3 – 14.6"	10.6 – 15"	10.6 – 13.8"	13.4 – 20.9"	60"	90° - 140°

†Upholstery colors	Red/ Grey	Blue/ Grey	Green/ Grey	Magenta/ Grey	Grey/ Grey
Size 1	£ 6242	£ 6243	£ 6244	£ 6241	£ 6240
Size 2	£ 6442	£ 6443	£ 6444	£ 6441	£ 6440

†Head rests size 1&2 (Required option) HCPCS E0955

£ 62611 - Shell \$255

£ 62621 - Neck support \$255

£ 6267 - Adjustable laterals \$305

£ 6268 - fitting, ball joint \$79



## Step 2a – Outdoor Frame Option (Required options are indicated below by † symbol)

**Advantage-chassis** – Folding aluminum frame, tilt-in-space, back frame suspension, face-to-face option, 0° safety lock out, cargo basket, unfold lock



	HCPCS	SIZE 1	SIZE 2
Advantage-chassis	E1234	£ 62023 \$1975	£ 64023 \$2025
†Castor pin lock, pair	K0073	67830 \$92	67830 \$92
†Bus tie-down loop kit (*)	K0108	67839 \$175	67839 \$175
Pneumatic wheels, full set	K0070	£ 66461 \$205	£ 66461 \$205
Frame padding, reversible	K0108	£ 6622 \$105	£ 6622 \$105
Umbrella w/ Adv.-chassis	K0108	£ 6594 \$10	£ 6594 \$10
Tip assist	K0108	£ 67502 \$69	£ 67502 \$69
Cargo basket	K0108	NC	NC

EASyS Modular S Technical Specs	SIZE 1	SIZE 2
Seat capacity	77.2 lbs	88.2 lbs
Frame capacity	110.2 lbs	121.3 lbs
Unit weight	41 lbs	48 lbs
Max seat tilt	35°	
Frame suspension	Yes	
Upholstery	Trevira-CS	
Foot angle	78° - 100°	
Size of foot rest	11" x 7.1"	12.6" x 8.3"
Floor to handle	28.2 – 46.1"	28.2 – 46.1"
Floor to seat	19"	19.5"
Max width	23.4"	27.2"
Fold length	33.5"	34"
Fold height	23.6"	27.2"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7.1" solid	
Rear wheels	10" solid	



\* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement. ExoMotion LLC dba Thomashilfen | [www.thomashilfen.us](http://www.thomashilfen.us) Effective Date January 1<sup>st</sup> 2018

**Step 2b – Indoor Frame Option (Required options are indicated below by † symbol)**

**Q-chassis** - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas springs, lockable rear wheels



	HCPCS	SIZE 1 & 2	
Q-chassis indoor hi/low	E1399	£ 6666	\$1399
†EASyS Adapter for Q-chassis	K0108	£ 62951	\$95
Push bar for Q-chassis	K0108	£ 6660	\$275

**Q-chassis Technical Specs**

Capacity	88.2 lbs
Seat tilt	-7 to +31°
Floor to seat	11.4 - 27.3"
Seat direction	front facing

**Step 3 – Positioning & General Options (Required options are indicated below by † symbol)**

Lower Body Positioning	HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle	E0952	Pair		£ 6648		\$30
Ankle Huggers, Bodypoint®	K0108	Pair		£ 213151		\$205
Footrest pad, height 1¾"	K0108	Each	£ 10622		£ 10023	\$95

Mid Body Positioning	HCPCS	Unit	Small	Medium	Large	Price
Abduction block, padded	E0957	Each	£ 6620		£ 6720	\$138
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420	£ 8470	\$71
	4 point, padded, SR buckle	E0978	Each	£ 8421	£ 8471	\$170
	Soft pelvic harness	E0978	Each	£ 6614	£ 6714	\$139
Crotch strap (fits 6629,8420,8470,8421,8471)	E0978	Each		£ 6341		\$25
Tray, height adjustable, clear with picture window	E0950	Each	£ 66281		£ 67281	\$190
Mounting hardware tray – †required when ordering 6628/6728	E1028	Each	£ 66282		£ 67282	\$205
Arm rest, height, depth and angle adjustable	E0973	Pair	£ 8431		£ 8478	\$375

Upper Body Positioning	HCPCS	Unit	Small	Medium	Large	Price	
Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385		£ 6585	\$145
	Vest style	E0960	Each	£ 6631	£ 6632	£ 6731	\$145
	H-style, posture support, AEL	E0960	Each		£ 14813		\$185
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		£ 214270		\$185
	PivotFit shoulder harness, Bodypoint®	E0960	Each		£ 214250		\$165
Trunk supports	Rigid lateral	E0956	Pair	£ 6634		£ 6734	\$200
	Flexible lateral w/ chest strap	E0956	Pair		£ 6633		\$285

General Accessories	HCPCS	Unit	Small	Medium	Large	Price	
Grab bar, padded, black	K0108	Each	£ 6621		£ 6721	\$105	
Footrest tube upholstery, grey	K0108	Pair		£ 6230		\$16	
Vent tray	E1029	Each	£ 6294		£ 6494	\$410	
Medical necessities pack, black	K0108	Each		£ 6639		\$67	
Oxygen tank holder, E-tank	E2208	Each		£ 665012		\$140	
Straps	Secure folded seat/frame	K0108	Each		£ 667011	\$25	
	Wrist "runaway" tether	K0108	Each		£ 667111	\$25	
Weather protection	Canopy, sun/rain with window, black	K0108	Each		£ 663761	\$295	
	Lower body windbreaker, black for 663761	K0108	Each		£ 661852	\$95	
	Canopy, rain only, clear	K0108	Each		£ 6638	\$195	
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		£ 6739	\$109	
	Rain cape, grey	K0108	Each	£ 6340		£ 6540	\$149
	Sleeping bag, summer weight, black	K0108	Each	£ 6619		£ 6719	\$295
	Sleeping bag, winter weight, black	K0108	Each	£ 6618		£ 6718	\$299
Umbrella, angle adjustable	K0108	Each		£ 6594		\$85	

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.