

EASyS Modular S – 2022 US Retail prices

info@thomashilfen.com | Phone 866.870.2122 | Fax 866.870.0801

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat Sizes (Required options are indicated below by † symbol)

£ Modular S seat £ Modular S seat



	size 1	size 2	HCPCS	Description		
†Seat	661011	\$420	671011	\$420	E2292	Solid contoured, integrated seat/back, 3 different head rest options for more head control; 140° back recline, growth adjustable; seat depth/ width adjustable from underneath, articulated foot plate height and angle adjustable; front/rear facing seat, easy to change direction; ventilated/washable upholstery ; knee angle adjustment with quick release lever; cozy comfortable extra seat cushion; calf panel; 5 different color options; tilt-in-space chassis available
†Back	661012	\$518	671012	\$518	E2293	
†Recline	661013	\$415	671013	\$415	E1014	
†Hip sup.	661014	\$230	671014	\$230	E0956	
†Head rest bracket	661015	\$170	671015	\$170	E1028	
†Calf panel	6635	\$45	6735	\$45	E0995	
†Foot plate, flip up	661016	\$90	671016	\$90	K0040	
†Legrest, elevating	661017	\$200	671017	\$200	E0990	

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Back Height#	Max Height	Hip Angle
Size 1	7.3 – 12"	7.5 – 11.8"	7.5 – 10.2"	11 – 16.5"	50"	90° - 140°
Size 2	7.3 – 14.6"	10.6 – 15"	10.6 – 13.8"	13.4 – 20.9"	60"	90° - 140°

†Upholstery colors	Red/ Grey	Blue/ Grey	Green/ Grey	Magenta/ Grey	Grey/ Grey
Size 1	£ 6242	£ 6243	£ 6244	£ 6241	£ 6240
Size 2	£ 6442	£ 6443	£ 6444	£ 6441	£ 6440

†Head rests size 1&2 (Required option) HCPCS E0955

£ 62611 - Shell \$285

£ 62621 - Neck support \$285

£ 6267 - Angle adjust. laterals \$335 £ 6268 - fitting, ball joint \$90



Step 2a – Outdoor Frame Option (Required options are indicated below by † symbol)

Advantage-chassis - Folding aluminum frame, tilt-in-space, back frame suspension, face-to-face option, 0° safety lock out, cargo basket, unfold lock



	HCPCS	SIZE 1	SIZE 2
Advantage-chassis	E1234	£ 62023 \$2295	£ 64023 \$2295
†Castor pin lock, pair	K0073	67830 \$96	67830 \$96
†Bus tie-down loop kit (*)	K0108	67839 \$205	67839 \$205
Pneumatic wheels, full set	K0070	£ 66461 \$249	£ 66461 \$249
Frame padding, black	K0108	£ 6622 \$115	£ 6622 \$115
Frame padding, red/ black	K0108	£ 66221 \$115	£ 66221 \$115
Tip assist/ curb tippers	K0108	£ 67502 \$89	£ 67502 \$89
Under seat storage basket	K0108	NC	NC



EASyS Modular S Technical Specs	SIZE 1	SIZE 2
Seat capacity	77.2 lbs	88.2 lbs
Frame capacity	110.2 lbs	121.3 lbs
Unit weight	41 lbs	48 lbs
Max seat tilt	35°	
Frame suspension	Yes	
Upholstery	Trevira-CS	
Foot angle	78° - 100°	
Size of foot rest	11" x 7.1"	12.6" x 8.3"
Floor to handle	28.2 – 46.1"	28.2 – 46.1"
Floor to seat	19"	19.5"
Max width	23.4"	27.2"
Fold length	33.5"	34"
Fold height	23.6"	27.2"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7.1" solid	
Rear wheels	10" solid	

* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement. ExoMotion LLC dba Thomashilfen | www.thomashilfen.us Effective Date January 1st 2022

Step 2b – Indoor Frame Option (Required options are indicated below by † symbol)

Q-chassis - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas springs, lockable rear wheels



	HCPCS	SIZE 1 & 2	
Q-chassis indoor hi/low, pedal controlled	E1399	£ 6666	\$1999
†EASyS Adapter for Q-chassis	K0108	£ 62951	\$130
Push bar for Q-chassis	K0108	£ 6660	\$298

Q-chassis Technical Specs

Capacity	88.2 lbs
Seat tilt	-7 to +31°
Floor to seat	11.4 - 27.3"
Seat direction	front facing

Step 3 – Positioning & General Options (Required options are indicated below by † symbol)

Lower Body Positioning	HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle	E0952	Pair		£ 6648		\$35
Ankle Huggers, Bodypoint®	K0108	Pair		£ 213151		\$205
Footrest pad, height 1¾"	K0108	Each	£ 10622		£ 10023	\$99

Mid Body Positioning						
	HCPCS	Unit	Small	Medium	Large	Price
Abduction block, padded	E0957	Each	£ 6620		£ 6720	\$160
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420	£ 8470	\$95
	4 point, padded, SR buckle	E0978	Each	£ 8421	£ 8471	\$135
	Soft pelvic harness/ leg harness	E0978	Each	£ 6614	£ 6714	\$175
Crotch strap (fits 6629,8420,8470,8421,8471)	E0978	Each		£ 6341		\$33
Tray, height adjustable, clear with picture window	E0950	Each	£ 66281		£ 67281	\$190
Mounting hardware tray – †required when ordering 6628/6728	E1028	Each	£ 66282		£ 67282	\$225

Upper Body Positioning							
	HCPCS	Unit	Small	Medium	Large	Price	
Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385		£ 6585	\$189
	Vest style	E0960	Each	£ 6631	£ 6632	£ 6731	\$189
	H-style, posture support, small, Bodypoint®	E0960	Each		£ 214806		\$199
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		£ 214270		\$199
	PivotFit shoulder harness, Bodypoint®	E0960	Each		£ 214250		\$189
Trunk supports	Rigid lateral, left & right, height adj.	E0956	Pair	£ 6634		£ 6734	\$245
	Flexible lateral l & r, w/chest strap, height adj.	E0956	Pair		£ 6633		\$295

General Accessories							
	HCPCS	Unit	Small	Medium	Large	Price	
Grab bar, padded, black	K0108	Each	£ 6621		£ 6721	\$125	
Footrest tube upholstery, grey	K0108	Pair		£ 6230		\$28	
Vent tray	E1029	Each	£ 6294		£ 6494	\$450	
Medical necessities pack, black	K0108	Each		£ 6639		\$80	
Oxygen tank holder, E-tank	E2208	Each		£ 665012		\$210	
Straps	Secure folded seat/frame	K0108	Each		£ 667011	\$30	
	Wrist "runaway" tether	K0108	Each		£ 667111	\$30	
Weather protection	Canopy, sun/rain with window, black	K0108	Each		£ 663761	\$325	
	Lower body windbreaker, black for 663761	K0108	Each		£ 661852	\$99	
	Canopy, rain only, clear	K0108	Each		£ 6638	\$245	
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		£ 6739	\$125	
	Rain cape, grey	K0108	Each	£ 6340		£ 6540	\$163
	Sleeping bag, summer weight, black	K0108	Each	£ 6619		£ 6719	\$299
	Sleeping bag, winter weight, black	K0108	Each	£ 6618		£ 6718	\$299
Umbrella, angle adjustable	K0108	Each		£ 6594		\$95	

*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.