


# EASyS Modular – 2018 US Retail prices

info@thomashilfen.com  
Phone 866.870.2122 | Fax 866.870.0801



Date	Buyer	PO#	Mark For
Address			
Contact	Phone ( )	Fax ( )	

## Step 1 – Seat Sizes (Required options are indicated below by † symbol)



	EASyS Modular seat size		EASyS Modular seat size 2		HCPCS	Description
†Seat	661001	\$475	671001	\$525	E2294	Solid contoured, integrated seat/back, 140° recline, growth adjustable; articulated footrest height and angle adjustable; front/rear facing seat; height and angle adjustable padded arm supports; height, depth and angle adjustable swing-away head rest bracket; 3 different head rest shapes; ventilated/ washable upholstery
†Back	661002	\$475	671002	\$525	E2293	
†Recline	661003	\$250	671003	\$250	E1014	
†Hip sup.	661004	\$200	671004	\$200	E0956	
†Bracket	661005	\$125	671005	\$125	E1028	
†Calf strap	6635	\$35	6735	\$35	E0995	

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Back height#	Max Height	Hip Angle
Size 1	6¾ – 11¾"	6¼ – 11"	7⅞ – 11¾"	10¼ – 16½"	50"	90° - 140°
Size 2	9¾ – 14¼"	9¾ – 13¾"	10¼ – 13¾"	12⅝ – 20½"	60"	90° - 140°

†Trunk pad colors	Purple	Classic	Berry	Pink	†Head rest	\$245	HCPCS E0955
Size 1	£ 6265	£ 6264	£ 6263	£ 6266		6260/6460	6261
Size 2	£ 6465	£ 6464	£ 6463	£ 6466		£ Standard	£ Neck supp.

## Step 2 – Frame Options (Required options are indicated below by † symbol)

**A-chassis** - Folding aluminum frame, tilt-in-space, dual spring suspension, cargo basket, tilt-in-space lever at push bar

**Q-chassis** - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas spring, lockable rear wheels



	HCPCS	SIZE 1	SIZE 2
A-chassis	E1234	£ 6202 \$1975	£ 6402 \$2025
†Castor pin lock, pair	K0073	67830 \$92	67830 \$92
†Bus tie-down loop kit (*)	K0108	67839 \$175	67839 \$175
Pneumatic wheels, full set	K0070	£ 6646 \$205	£ 6646 \$205
Frame padding, reversible	K0108	£ 6622 \$105	£ 6622 \$105
Umbrella w/ A-chassis	K0108	£ 6594 \$10	£ 6594 \$10
Tip Assist	K0108	£ 6750 \$69	£ 6750 \$69
Cargo basket	K0108	NC	NC

	HCPCS	SIZE 1 & 2
Q-chassis	E1399	£ 6666 \$1399
Push bar for Q-chassis	K0108	£ 6660 \$275

### Q-chassis Technical Specs \*\*

Capacity	88 lbs
Seat tilt	-7 to +31°
Floor to seat	11 - 28"
Seat direction	front facing

### EASyS Modular Technical Specs \*\*

	SIZE 1	SIZE 2
Seat capacity	77 lbs	88 lbs
Frame capacity	99 lbs	110 lbs
Unit weight	36 lbs	42 lbs
Max seat tilt	35°	
Seat suspension, lockable	Yes	
Upholstery	Trevira-CS	
Height adjustment arm support	3¾ – 5½"	4¾ – 7"
Foot angle	78° - 100°	
Size of foot rest	11½" x 7⅞"	13" x 8¾"
Floor to handle	27 – 46"	31 – 48"
Floor to seat	19"	19½"
Max width	23"	27"
Fold length	32"	35½"
Fold height	18½"	18½"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7½" solid	
Rear wheels	10" solid	

**Z-chassis** - indoor base/feeder/therapy frame, manual height/ tilt adjustable, lockable rear wheels



	HCPCS	SIZE 1 & 2
Z-chassis	E1399	£ 6691 \$899

### Z-chassis Technical Specs \*\*

Capacity	88 lbs
Seat tilt	-15 to +35°
Floor to seat	15 - 23"
Seat direction	front facing

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

**Step 3 – Positioning & General Options (Required options are indicated below by † symbol)**

<b>Lower Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle		E0952	Pair		£ 6648		\$30
Ankle Huggers, Bodypoint®		K0108	Pair		£ 213151		\$205
Footrest pad, height 1¾"		K0108	Each	£ 10622		£ 10023	\$95

<b>Mid Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Abduction block, padded		E0957	Each	£ 6620		£ 6720	\$138
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420		£ 8470	\$71
	4 point, padded, SR buckle	E0978	Each	£ 8421		£ 8471	\$170
	Soft pelvic harness	E0978	Each	£ 6614		£ 6714	\$139
Crotch strap (fits 6629,8420,8470,8421,8471)		E0978	Each		£ 6341		\$25
Tray, height adjustable, clear with picture window		E0950	Each	£ 66281		£ 67281	\$190
Mounting hardware tray – †required when ordering 6628/6728		E1028	Each	£ 66282		£ 67282	\$205

<b>Upper Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385		£ 6585	\$145
	Vest style	E0960	Each	£ 6631	£ 6632	£ 6731	\$145
	H-style, posture support, AEL	E0960	Each		£ 14813		\$185
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		£ 214270		\$185
	PivotFit shoulder harness, Bodypoint®	E0960	Each		£ 214250		\$165
Trunk supports	Rigid lateral	E0956	Pair	£ 6634		£ 6734	\$200
	Flexible lateral w/ chest strap	E0956	Pair		£ 6633		\$285

<b>General Accessories</b>		HCPCS	Unit	Small	Medium	Large	Price
Grab bar, padded, black		K0108	Each	£ 6621		£ 6721	\$105
Vent tray		E1029	Each		£ 110004		\$410
Oxygen tank holder, E-tank		E2208	Each		£ 665012		\$140
Straps	Secure folded seat/frame	K0108	Each		£ 667011		\$25
	Wrist "runaway" tether	K0108	Each		£ 667111		\$25
Weather protection	Canopy, sun/rain with window, black	K0108	Each		£ 663761		\$295
	Lower body windbreaker, black for 663761	K0108	Each		£ 661852		\$95
	Canopy, rain only, clear	K0108	Each		£ 6638		\$195
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		£ 6739		\$109
	Rain cape, grey	K0108	Each	£ 6340		£ 6540	\$149
	Sleeping bag, summer weight, black	K0108	Each	£ 6619		£ 6719	\$295
	Sleeping bag, winter weight, black	K0108	Each	£ 6618		£ 6718	\$299
Umbrella, angle adjustable		K0108	Each		£ 6594		\$85

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

# Back height is without headrest  
 \*\*More details in Technical Manual