

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat & Frame (Required options are indicated below by † symbol)



	£ Swifty	HCPCS	Description
†Seat	680621 \$1995	E1236	Folding frame with integrated/ contoured seat, fixed tilt, adjustable seat depth from underneath, footrest height adjustable, back recline, suitable for Early Intervention use
†Back	680622 \$90	E1014	
†Castor locks	6810 \$96	K0073	Castor pin lock, pair
†Bus tie-down	6830 \$205	K0108	Bus tie-down loop kit *
H-belt	6629 \$30	E0960	With purchase of Swiftly

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Seat to Head	Max Height	Hip Angle
Swiftly	6¼" – 13"	8½" – 11¼"	13½"	16" – 24½"	40¼"	90° to 115°

Step 2 – Positioning & General Options

Lower Body Positioning	HCPCS	Unit	Small >>> Medium >>> Large	Price
Foot straps with cam buckle	E0952	Pair	£ 6648	\$35
Ankle Huggers**	K0108	Pair	£ 213151	\$205
Footrest pad, height 1¾"	K0108	Each	£ 10023	\$99

Mid Body Positioning

Abduction block, padded	E0957	Each	£ 6620	\$160	
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420	\$95
	Soft pelvic harness	E0978	Each	£ 6614 £ 6714	\$158
Crotch strap (fits 6629,8420)	E0978	Each	£ 6341	\$33	
Tray, height adjustable, clear with picture window	E0950	Each	£ 68221	\$190	
Mounting hardware tray †required ordering 6822	E1028	Each	£ 68222	\$225	

Upper Body Positioning

Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385 £ 6585	\$189
	Vest style	E0960	Each	£ 6631 £ 6632	\$189
	H-style, w. integrated hip belt	E0960	Each	£ 66297	\$110
	H-style, posture support, small**	E0960	Each	£ 214806	\$199
	PivotFit shoulder harness, dynamic**	E0960	Each	£ 214270	\$199
	PivotFit shoulder harness**	E0960	Each	£ 214250	\$189
Trunk supports	Rigid lateral, padded, l&r, height adj.	E0956	Pair	£ 6634	\$245
	Flexible lateral with chest strap	E0956	Pair	£ 6633	\$295
Head supports	Contoured pillow	E0955	Each	£ 6824	\$175

Swiftly	Technical Specs #
Capacity	77 lbs
Unit weight	27 lbs
Seat tilt	15°
Back recline	90-115°
Upholstery	Trevira-CS
Floor to handle	28"-45"
Floor to seat	19"
Max width	24"
Fold length	28¾"
Fold height	15"
Brakes	foot operated
Seat direction	front facing
Front wheels	7½" solid
Rear wheels	10" solid

More details in Technical Manual



General Accessories

Seat minimizer kit for Early Intervention use	K0108	Set	£ 6834	\$290	
Pneumatic wheels, full set	K0070	Set	£ 6646	\$245	
Under seat storage basket	K0108	Each	£ 6813	\$75	
Frame padding, black	K0108	Pair	£ 6622	\$115	
Frame padding, red/ black, reversible	K0108	Pair	£ 66221	\$115	
Grab bar, padded, black	K0108	Each	£ 6833	\$125	
Tip assist	K0108	Each	£ 6829	\$89	
Straps	Secure folded seat/frame	K0108	Each	£ 667011	\$30
	Wrist "runaway" tether	K0108	Each	£ 667111	\$30
Weather protection	Canopy, sun/rain with window, black	K0108	Each	£ 663761	\$325
	Lower body windbreaker for 663761	K0108	Each	£ 661852	\$99
	Canopy, rain only, clear	K0108	Each	£ 6638	\$245
	Rain shield, clear for canopy 663761	K0108	Each	£ 6739	\$125
	Rain cape, grey	K0108	Each	£ 6340	\$163
	Sleeping bag, summer weight, black	K0108	Each	£ 6619	\$299
	Sleeping bag, winter weight, black	K0108	Each	£ 6618	\$299
	Umbrella, angle adjustable, grey	K0108	Each	£ 6594	\$95

* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle. ** Bodypoint®