

# Swiftly – Rehab stroller

All details online:  
[www.thomashilfen.us/swifty1](http://www.thomashilfen.us/swifty1)  or [.../swifty2](http://www.thomashilfen.us/swifty2) 

## Accessories (selection)



Sleeping bag, padded (summer) or woven fur (winter) (shown on size 1)



Sun and rain cover with leg blanket (shown on size 1)



Umbrella (shown on size 1)




Basket (shown on size 1)



Butterfly chest harness (flexible), abduction block and grip rail (shown on size 2)



Head pillow, lateral trunk support with chest belt and 2-point pelvic harness (shown on size 2)



Soft pelvic harness (shown on size 2)



Tip assist (shown on size 2)

You will find the total overview of all rehab accessories on our website: [www.thomashilfen.us](http://www.thomashilfen.us)

Technical data	Swiftly 1	Swiftly 2
Seat depth	22 (17 cm) - 28.5 cm / (with seat minimizer) 8.7 (6.9) - 11.2"	34 cm (24.5 cm) - 40 cm / 13.4" (9.6") - 15.7"
Seat width	34 cm (23 cm) / (with seat minimizer) 13.4" (9.1")	39.5 cm (29.5 cm) / 15.6" (11.6")
Back height	62 cm / 24.4"	70 cm / 27.6"
min. shoulder height	37 cm (26 cm) / (with seat minimizer) 14.6" (10.2")	47 cm (30 cm) / 18.5" (11.8")
Lower leg length	16 (12 cm) - 33 cm / (with seat minimizer) 6.3 (4.7") - 13"	21 (21 cm) - 37 cm / 8.3 (8.3 ") - 14.6"
Hip angle	90° to 115°	90° to 100°
Foot angle	90°	90°
Seat tilt-in-space fix	+15°	+20°
Size of foot rest (w x d)	31 x 18 cm / 12.2 x 7.1"	31 x 20 cm / 12.2 x 7.9"

Technical data	Swiftly 1	Swiftly 2
Push bar height	80 - 119 cm / 31.5 - 46.9"	82 - 124 cm / 32.3 - 48.8"
Wheel size swivel wheels (front)	18 cm / 7.1"	18 cm / 7.1"
Wheel size (rear)	25 cm / 9.8"	25 cm / 9.8"
Overall dimensions (l x w x h)	61 x 102 x 98 cm / 24 x 40.2 x 38.6"	67 x 113 x 98 cm / 26.4 x 44.5 x 38.6"
Folded dimensions (l x w x h)	61 x 38 x 73 cm / 24 x 15 x 28.7"	67 x 49 x 79 cm / 26.4 x 19.3 x 31.1"
Weight	12.4 kg / 27.3 lb	18.6 kg / 41 lb
max. user weight	35 kg / 77.2 lb	50 kg / 110.2 lb



**NEW** FOR SIZE 2:  
 Seat depth adjustable without tools

# Swiftly

## The lightweight rehab strollers for both small & big kids



Thomashilfen North America (ExoMotion, LLC) · 14221 29th St E, #101 · Sumner, WA 98390 · USA  
 Phone: 866.870.2122 (toll free) · Fax: 866.870.0801 (toll free) · info@thomashilfen.com · www.thomashilfen.us



Item code: 9800 95773 0000, status: 2023-03

# Simple use, light & swift: Here comes **Swiftly** ...

“ A stable, comfortable rehab stroller with a fixed seat; it can withstand the many stresses of day-to-day life and is quickly ready for use. Now available in two sizes for growing children. ”

Comfortable sitting for the child: the washable seat cover is soft and breathable – preventing perspiration and providing a thoroughly cozy feeling



Back angle can be adjusted for positioning and relaxed seating



Fold lock (only for size 2) keeps everything securely folded together.

**NEW** FOR SIZE 2:  
Seat depth adjustable without tools



 size 2

 size 1



With the spacious Swiftly size 2

- big kids that weigh up to 50 kg / 110.2 lb
- and are up to 150 cm / 59.1" tall can be cared for with a rehab stroller



Transit tie-down-kit for safe transport in a vehicle (available for both sizes, standard for size 2)

Small folding size (even size 2) – Especially quick and easy to fold with a flick of the wrist.



Foldable and height-adjustable footrest – Children can get in and out of the stroller independently.



Seat depth adjustable without tools – even while the child is sitting in the rehab stroller  
**NOW available also for SIZE 2**



Seat, back and footrest provide safe and firm support – a convenient place to sit for the child



H-belt with padded harness as a standard feature



Seat minimizer for smaller children (accessory) – grows with your child for a very long time

## Swiftly & Swiftly 2

The lightweight rehab strollers for both small & big kids