

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat & Frame (Required options are indicated below by † symbol)



	£ Swiftly 2	HCPCS	Description
†Seat	680821 \$2825	E1236	Folding frame with integrated/ contoured seat, fixed tilt, adjustable seat depth from underneath, footrest height adjustable, back recline, suitable for Early Intervention use
†Back	680822 \$125	E1014	
†Castor locks	6810 \$96	K0108	Castor pin locks, pair
†Bus tie-down	6830 \$175	K0108	Bus tie-down loop kit *
H-belt	6629 \$5	E0960	With purchase of Swiftly
Seat Size Chart			
	Lower leg length	Seat depth	Seat width
Swiftly 2	8.3" – 14.6"	14.8"	15.6"
	Back height	Min. shoulder height	Hip angle
Swiftly 2	27.6"	11.8"	90° to 100°

Step 2 – Positioning & General Options

Lower Body Positioning	HCPCS	Unit	Price
Foot straps with cam buckle	E0952	Pair	£ 6648 \$30
Ankle Huggers**	K0108	Pair	£ 213151 \$205
Footrest pad, height 1¾"	K0108	Each	£ 10023 \$95

Mid Body Positioning	HCPCS	Unit	Price
Abduction block, padded	E0957	Each	£ 6720 \$138
Hip supports	2-point, padded, SR buckle	E0978	Each £ 8470 \$80
	4-point, padded, SR buckle	E0978	Each £ 8471 \$135
	Soft pelvic harness	E0978	Each £ 6714 \$145
Crotch strap (fits 6629,8470)	E0978	Each	£ 6341 \$28
Tray, height adjustable, clear with picture window	E0950	Each	£ 67281 \$190
Mounting hardware tray †required ordering 6728	E1028	Each	£ 67282 \$205

Upper Body Positioning	HCPCS	Unit	Price
Chest harness supports	Butterfly style, neoprene	E0960	Each £ 6585 \$165
	Vest style	E0960	Each £ 6632 \$165
	H-style, w. integrated hip belt	E0960	Each £ 66297 \$96
	H-style, posture support, small**	E0960	Each £ 214806 \$199
	PivotFit shoulder harness, dynamic**	E0960	Each £ 214270 \$199
	PivotFit shoulder harness**	E0960	Each £ 214250 \$189
Trunk supports	Rigid lateral	E0956	Pair £ 6734 \$230
	Flexible lateral with chest strap	E0956	Pair £ 6633 \$285
Head supports	Contoured pillow, black	E0955	Each £ 6837 \$140
	Contoured pillow, red	E0955	Each £ 68371 \$140

General Accessories			
Seat minimizer kit for Early Intervention use	K0108	Set	£ 6836 \$275
Seat width reducer, black (seat width 11.6")	K0108	Pair	£ 68410 \$130
Seat width reducer, red (seat width 11.6")	K0108	Pair	£ 6841 \$130
Pneumatic wheels, full set	K0070	Set	£ 6646 \$215
Under seat storage basket	K0108	Each	£ 6835 \$70
Frame padding, black	K0108	Pair	£ 6622 \$105
Frame padding, red/black, reversible	K0108	Pair	£ 66221 \$105
Grab bar, padded, black	K0108	Each	£ 6721 \$110
Tip assist	K0108	Each	£ 6829 \$75
Straps	Secure folded seat/frame	K0108	Each £ 667011 \$25
	Wrist "runaway" tether	K0108	Each £ 667111 \$25
Weather protection	Canopy, sun/rain with window, black	K0108	Each £ 663761 \$295
	Lower body windbreaker for 663761	K0108	Each £ 661852 \$95
	Canopy, rain only, clear	K0108	Each £ 6638 \$215
	Rain shield, clear for canopy 663761	K0108	Each £ 6739 \$109
	Rain cape, grey	K0108	Each £ 6540 \$149
	Sleeping bag, summer weight, black	K0108	Each £ 6719 \$295
	Sleeping bag, winter weight, black	K0108	Each £ 6718 \$299
Umbrella, angle adjustable, grey	K0108	Each £ 6594 \$85	

Swiftly 2	Technical Specs #
Capacity	110 lbs
Unit weight	40 lbs
Seat tilt	20°
Back recline	90-100°
Upholstery	Trevira-CS
Floor to handle	32.3"-48.8"
Floor to seat	
Max width	26.4"
Fold length	31.1"
Fold height	19.3"
Brakes	foot operated
Seat direction	front facing
Front wheels	7.1" solid
Rear wheels	9.8" solid



* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

** Bodypoint®