

# Swifty & Swifty 2

## TECHNICAL MANUAL



### Rehab Stroller

- Basic models
- Accessories



Swiftly



Swiftly 2

## Swiftly – The lightweight rehab stroller for children

You can also find more information on the Internet: [www.thomashilfen.com](http://www.thomashilfen.com)

| <b>Content</b>                   | <b>Page</b>  |
|----------------------------------|--|
| BASIC MODEL                      |  |
|                                  | Swiftly . . . . . 4  |
|                                  | Swiftly 2 . . . . . 5  |
| ACCESSORIES – Feet / Legs        |  |
|                                  | Foot straps . . . . . 6                                      |
| ACCESSORIES – Seat               |  |
|                                  | Abduction block . . . . . 6                                  |
|                                  | Grip rail . . . . . 6  |
|                                  | Tray . . . . . 6   |
|                                  | Lateral trunk support, rigid . . . . . 7                     |
|                                  | Lateral trunk support with chest belt . . . . . 7            |
|                                  | Lateral seat upholstery . . . . . 7                          |
|                                  | Frame padding . . . . . 7                                    |
| ACCESSORIES – Seat minimizer     |  |
|                                  | Seat minimizer – set . . . . . 8                             |
|                                  | Footpad (Seat minimizer) . . . . . 8                         |
|                                  | Pelvic upholstery (Sitzverkleinerer) . . . . . 9             |
|                                  | Seat width reducer (Seat minimizer) . . . . . 9              |
|                                  | Backrest upholstery (Seat minimizer) . . . . . 9             |
|                                  | Head pillow (Seat minimizer) . . . . . 9                     |
| ACCESSORIES – Belts              |  |
|                                  | 2-point pelvic harness with front-pull . . . . . 10          |
|                                  | 4-point pelvic harness with front-pull . . . . . 10          |
|                                  | Crotch strap . . . . . 10                                    |
|                                  | Soft pelvic harness . . . . . 10                             |
|                                  | Butterfly/chest harness . . . . . 11                         |
|                                  | Butterfly/chest harness, flexible . . . . . 11               |
|                                  | Seat waistcoat . . . . . 11                                  |
| ACCESSORIES – Head               |  |
|                                  | Head pillow . . . . . 11                                     |
| ACCESSORIES – Weather protection |  |
|                                  | Sun and rain cover with leg blanket . . . . . 12             |
|                                  | Sun and rain cover . . . . . 12                              |
|                                  | Rain shield for sun and rain cover, transparent . . . . . 12 |
|                                  | Rain cover, transparent, with leg blanket . . . . . 13       |
|                                  | Umbrella . . . . . 13  |
|                                  | Sleeping bag (padded + woven fur) . . . . . 13               |
| ACCESSORIES – Chassis            |  |
|                                  | Basket . . . . . 14  |
|                                  | Pneumatic wheels . . . . . 14                                |
|                                  | Tip assist . . . . . 15                                      |

# BASIC MODEL



## Swift

| Item code  | 9800068065000                                |
|--|--|
| Seat depth (with seat minimizer (accessory))       | 22 (17 cm) - 28,5 cm /<br>8.7 (6.7") - 11.2" |
| Seat width (with seat minimizer (accessory))       | 34 cm (23 cm) /<br>13.4" (9.1")              |
| Back height  | 62 cm /<br>24.4"                             |
| min. shoulder height (with seat minimizer)         | 37 cm (26 cm) /<br>14.6" (10.2")             |
| Lower leg length (with seat minimizer (accessory)) | 16 (12 cm) - 33 cm /<br>6.3 (4.7") - 13"     |
| Back angle   | 90° to 115°                                  |
| Foot angle   | 90°  |
| Seat tilt fix                                      | +15°   |
| Footrest size (W x L)                              | 31 x 18 cm /<br>12.2 x 7.1"                  |
| Overall dimensions (W x H x L)                     | 61 x 102 x 98 cm /<br>24.0 x 40.2 x 38.6"    |
| Dimensions when folded (W x H x L)                 | 61 x 38 x 73 cm /<br>24.0 x 15.0 x 28.7"     |
| Push bar height                                    | 80 - 119 cm /<br>31.5 - 46.9"                |
| Wheel size, swivel wheels (front)                  | 7,5"   |
| Wheel size (rear)                                  | 10"  |
| max. carrying capacity basket (accessory)          | 3 kg / 6.6 lbs                               |
| Weight   | 12.4 kg / 27.3 lbs                           |
| max. user weight                                   | 35 kg / 77.2 lbs                             |

### Standard Equipment

foldable chassis with integrated seat unit; seat and back pads; seat side guide pads; H-belt; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels with wheel lock; noflat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; crash eyes; bracket for sun and rain cover

### Material Specifications

Frame: powder-coated aluminium

Padding: 100 % Polyester

The material used for the padding is "flame retardant" and was tested in accordance with DIN EN 1021-1.

The material and components underwent standard testing to determine the extent of their resistance to inflammation and were assessed as being very unlikely to catch fire.



## Swift 2

| Item code  | 9800068080010                                 |
|--|---|
| Seat depth (with seat minimizer (accessory))           | 34 - 40 cm (24,5 cm) /<br>13.4 - 15.7" (9.6") |
| Seat width (with seat minimizer (accessory))           | 39,5 cm (29,5 cm) /<br>15.6" (11.6")          |
| Back height  | 70 cm / 27.6"                                 |
| min. shoulder height (with seat minimizer (accessory)) | 47 cm (30 cm) /<br>18.5" (11.8")              |
| Lower leg length                                       | 21 - 37 cm /<br>8.3 - 14.6"                   |
| Side guide (W x H)                                     | 25 x 16 cm /<br>9.8 x 6.3"                    |
| Back angle   | 90° - 100°                                    |
| Foot angle   | 90°   |
| Seat angle, fixed                                      | + 20°   |
| Safe tilting up to a maximum of                        | 10°   |
| Size of footrest (W x D)                               | 31 x 20 cm /<br>12.2 x 7.9"                   |
| Overall dimensions (L x W x H)                         | 98 x 67 x 113 cm /<br>38.6 x 26.4 x 44.5"     |
| Dimensions when folded (L x W x H)                     | 79 x 67 x 49 cm /<br>31.1 x 26.4 x 19.3"      |
| Push bar height  | 82 - 124 mm /<br>32.3 - 48.8"                 |
| Turning circle   | 126 cm /<br>49.6"                             |
| Wheel size (front (swivelling wheels) / rear)          | 18 cm / 25 cm /<br>7.1" / 9.8"                |
| Wheel width (front (swivelling wheels) / rear)         | 4,5 cm / 5 cm /<br>1.8" / 2"                  |
| Weight   | 18,6 kg / 41 lbs                              |
| Max. load capacity of the basket (accessory)           | 3 kg / 6.6 lbs                                |
| Max. user weight                                       | 50 kg / 110.2 lbs                             |

### Standard Equipment

foldable chassis with integrated seat unit; seat and back pads; seat side guide pads; H-belt; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels with wheel lock; noflat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; crash eyes; bracket for sun and rain cover; folding lock

### Material Specifications

Frame: powder-coated aluminium

Padding: Polyester – Trevira CS

The material used for the padding is "flame retardant" and was tested in accordance with DIN EN 1021-1 and DIN EN 1021-2.


The material and components underwent standard testing to determine the extent of their resistance to inflammation and were assessed as being very unlikely to catch fire.

The 100% Trevira CS material is a high-quality, breathable material and tested positive for biocompatibility. (DIN EN ISO 10993-5 and 10993-10 Trevira CS, black)

# ACCESSORIES



| Foot straps        | Swiftly              | Swiftly 2            |
|--------------------|----------------------|----------------------|
| <b>Item code</b>   | <b>9800066482000</b> | <b>9800066482000</b> |
| <b>Belt length</b> | 87 cm / 34.3"        | 87 cm / 34.3"        |
| <b>Colour</b>      | black                | black                |



| Abduction block         | Swiftly   | Swiftly 2                            |
|-------------------------|---|--------------------------------------|
| <b>Item code</b>        | <b>9800066206000</b>                            | <b>9800067206000</b>                 |
| <b>Size (w x h x l)</b> | 6 x 9 x 10 cm<br>2.4 x 3.5 x 3.9"               | 7 x 10,5 x 11 cm<br>2.8 x 4.1 x 4.3" |
| <b>Colour</b>           | grey  | grey                                 |
| <b>Note</b>             | cannot be used in combination with crotch strap |                                      |



| Grip rail                      | Swiftly                                 | Swiftly 2                                      |
|--------------------------------|---|--|
| <b>Item code</b>               | <b>9800068331000</b>                    | <b>9800067215000</b>                           |
| <b>Height to seat surface</b>  | 14,5 / 14 / 16,5 cm<br>5.7 / 5.5 / 6.5" | 14 / 17 / 19 / 22 cm<br>5.5 / 6.7 / 7.5 / 8.7" |
| <b>Width (inner dimension)</b> | 37 cm / 14.6"                           | 38 cm / 14.9"                                  |
| <b>Colour</b>                  | grey                                    | grey   |



| Tray                          | Swiftly                   | Swiftly 2                     |
|-------------------------------|---------------------------|-------------------------------|
| <b>Item code</b>              | <b>9800068220000</b>      | <b>9800067281010</b>          |
| <b>Size (w x d)</b>           | 56 x 33,2 cm<br>22 x 13"  | 60 x 36,2 cm<br>23.6 x 14.3"  |
| <b>Height to seat surface</b> | 15 - 22 cm / 5.9 - 8.7"   | 15 - 22 cm / 5.9 - 8.7"       |
| <b>Edge height (inside)</b>   | 2,5 cm / 1"               | 2,5 cm / 1"                   |
| <b>Angle adjustment</b>       | approx. 27°               | approx. 27°                   |
| <b>Depth adjustment range</b> | 16 cm / 6.3"              | 19 cm / 7.5"                  |
| <b>Cut-out (w x d)</b>        | 41 x 14 cm<br>16.1 x 5.5" | 45,5 x 17,5 cm<br>17.9 x 6.9" |
| <b>Material thickness</b>     | 8 mm / 0.3"               | 8 mm / 0.3"                   |
| <b>Colour tray surface</b>    | transparent / colourless  | transparent / colourless      |



## Lateral trunk support, rigid

|                                   | Swiftly                 | Swiftly 2               |
|-----------------------------------|-------------------------|-------------------------|
| <b>Item code</b>                  | <b>9800066346000</b>    | <b>9800067346000</b>    |
| <b>Supporting surface (h x d)</b> | 15 x 13 cm / 5.9 x 5.1" | 17 x 15 cm / 6.7 x 5.9" |
| <b>Colour</b>                     | grey                    | grey                    |



## Lateral trunk support, flexible, with chest belt

|                                   | Swiftly                 | Swiftly 2               |
|-----------------------------------|-------------------------|-------------------------|
| <b>Item code</b>                  | <b>9800066332000</b>    | <b>9800066332000</b>    |
| <b>Supporting surface (h x d)</b> | 10 x 20 cm / 3.9 x 7.9" | 10 x 20 cm / 3.9 x 7.9" |
| <b>Colour</b>                     | grey                    | grey                    |



## Lateral seat upholstery

|                     | Swiftly                 | Swiftly 2                |
|---------------------|-------------------------|--------------------------|
| <b>Item code</b>    | <b>5005100228000</b>    | <b>9700442011000</b>     |
| <b>Size (h x d)</b> | 15 x 24 cm / 5.9 x 9.4" | 16 x 26 cm / 6.2 x 10.2" |
| <b>Colour</b>       | red                     | red                      |



## Frame padding

|                                 | Swiftly                                  | Swiftly 2                                |
|---------------------------------|--|--|
| <b>Item code</b>                | <b>9800066226000</b>                     | <b>9800066226000</b>                     |
| <b>Outer dimensions (w x h)</b> | (10 - 15) x 45 cm<br>(3.9 - 5.9) x 17.7" | (10 - 15) x 45 cm<br>(3.9 - 5.9) x 17.7" |
| <b>Colour</b>                   | black                                    | black                                    |

|                                 | Swiftly                                  | Swiftly 2                                |
|---------------------------------|--|--|
| <b>Item code</b>                | <b>9800066227000</b>                     | <b>9800066227000</b>                     |
| <b>Outer dimensions (w x h)</b> | (10 - 15) x 45 cm<br>(3.9 - 5.9) x 17.7" | (10 - 15) x 45 cm<br>(3.9 - 5.9) x 17.7" |
| <b>Colour</b>                   | red / black                              | red / black                              |



# ACCESSORIES



## Seat minimizer – set **Swifty** (consists of: head pillow, backrest upholstery, pelvic upholstery, footpad)

|   |                           |
|---|---------------------------|
| <b>Item code</b>                        | <b>9800068340000</b>      |
| <b>Head pillow height (bottom edge)</b> | 42 - 48 cm / 16.5 - 18.9" |
| <b>Thoraxunterstützung Rückenhöhe</b>   | 27 cm / 10.6"             |
| <b>Seat upholstery seat depth</b>       | 16 cm / 6.3"              |
| <b>Seat upholstery seat width</b>       | 22 cm / 8.6"              |
| <b>Seat upholstery height</b>           | 12 cm / 4.7"              |
| <b>Footpad height</b>                   | 4,5 cm / 1.8"             |
| <b>Colour</b>                           | black                     |

|                   |                      |
|-------------------|----------------------|
| <b>Item code</b>  | <b>9800068341000</b> |
| <b>Dimensions</b> | s. option black      |
| <b>Colour</b>     | red                  |

### Dimensional changes basic model with seat minimizer:

|                             |                            |
|-----------------------------|----------------------------|
| <b>Seat depth</b>           | 17 - 28,5 cm / 6.7 - 11.2" |
| <b>Seat width</b>           | 23 cm / 9.1"               |
| <b>min. shoulder height</b> | 26 cm / 10.2"              |
| <b>Lower leg length</b>     | 12 - 33 cm / 4.7 - 13"     |



## Seat minimizer set **Swifty 2** (consists of: head pillow, seat wedge, lateral seat guide padding (right and left))

|   |                              |
|---|------------------------------|
| <b>Item code</b>                                  | <b>9800068360000</b>         |
| <b>Head pillow height (bottom edge)</b>           | 48 - 57 cm /<br>18.9 - 22.4" |
| <b>Seat depth with seat wedge</b>                 | 24,5 cm / 9.6"               |
| <b>Seat width with lateral seat guide padding</b> | 29 cm / 11.4"                |
| <b>Lateral seat guide padding height</b>          | 16 cm / 6.2"                 |
| <b>Colour</b>                                     | black                        |

### Dimensional changes basic model with seat minimizer:

|                             |                               |
|-----------------------------|-------------------------------|
| <b>Seat depth</b>           | 24,5 - 30,5 cm /<br>9.6 - 12" |
| <b>Seat width</b>           | 29,5 cm / 11.6"               |
| <b>min. shoulder height</b> | 30 cm / 11.8"                 |



## Footpad **Swifty** –

(part of seat minimizer set)

|                     |                          |   |
|---------------------|--------------------------|---|
| <b>Item code</b>    | <b>5005100239005</b>     | – |
| <b>Size (w x d)</b> | 31 x 17 cm / 12.2 x 6.7" | – |
| <b>Height</b>       | 4,5 cm / 1.8 "           | – |
| <b>Colour</b>       | black                    | – |





| Pelvic upholstery<br>(part of seat minimizer set) | Swiftly | – |
|---|---------|---|
|---|---------|---|

|                                |                             |   |
|--------------------------------|-----------------------------|---|
| <b>Item code</b>               | <b>5005100229005</b>        | – |
| <b>Size upholstery (h x d)</b> | 11,5 x 21,5 cm / 4.5 x 8.5" | – |
| <b>Width upholstery</b>        | 4 cm / 1.6"                 | – |
| <b>Colour</b>                  | black                       | – |


|                   |                      |   |
|-------------------|----------------------|---|
| <b>Item code</b>  | <b>5005100227000</b> | – |
| <b>Dimensions</b> | s. option black      | – |
| <b>Colour</b>     | red                  | – |



| Seat width reducer<br>(part of seat minimizer set (lateral seat guide padding)) | – | Swiftly 2 |
|---|---|-----------|
|---|---|-----------|

|                                |   |                         |
|--------------------------------|---|-------------------------|
| <b>Item code</b>               | – | <b>9800068410000</b>    |
| <b>Size upholstery (h x d)</b> | – | 16 x 24 cm / 6.3 x 9.4" |
| <b>Width upholstery</b>        | – | 5 cm / 2"               |
| <b>Colour</b>                  | – | black                   |

|                   |   |                      |
|-------------------|---|----------------------|
| <b>Item code</b>  | – | <b>9800068411000</b> |
| <b>Dimensions</b> | – | s. option black      |
| <b>Colour</b>     | – | red                  |



| Backrest upholstery<br>(part of seat minimizer set) | Swiftly | – |
|---|---------|---|
|---|---------|---|

|                                    |                         |   |
|------------------------------------|-------------------------|---|
| <b>Item code</b>                   | <b>5005100219005</b>    | – |
| <b>Size (h x w)</b>                | 28 x 26 cm / 11 x 10.2" | – |
| <b>Depth</b>                       | 12,5 cm / 4.9"          | – |
| <b>Width</b>                       | 21,5 cm / 8.5"          | – |
| (between lateral trunk upholstery) |                         |   |
| <b>Colour</b>                      | black                   | – |

|                   |                      |   |
|-------------------|----------------------|---|
| <b>Item code</b>  | <b>5005100226000</b> | – |
| <b>Dimensions</b> | s. option black      | – |
| <b>Colour</b>     | red                  | – |



| Head pillow<br>(part of seat minimizer set) | Swiftly | Swiftly 2 |
|---|---------|-----------|
|---|---------|-----------|

|                         |                                       |   |
|-------------------------|---------------------------------------|---|
| <b>Item code</b>        | <b>5005100209005</b>                  | – |
| <b>Size (w x h x d)</b> | 26 x 13,5 x 7,5 cm<br>10.2 x 5.3 x 3" | – |
| <b>Colour</b>           | black                                 | – |

|                         |                      |                                      |
|-------------------------|----------------------|--------------------------------------|
| <b>Item code</b>        | <b>5005100225000</b> | <b>9800068420000</b>                 |
| <b>Size (w x h x d)</b> | s. option black      | 30,5 x 16 x 7,5 cm<br>12 x 6.3 x 3 " |
| <b>Colour</b>           | red                  | red                                  |


# ACCESSORIES




| 2-point pelvic harness with front-pull | Swiftly              | Swiftly 2            |
|--|----------------------|----------------------|
| <b>Item code</b>                       | <b>9800084201000</b> | <b>9800084701000</b> |
| <b>Upholstered area (w x h)</b>        | 15 x 5 cm / 5.9 x 2" | 17 x 5 cm / 6.7 x 2" |
| <b>Belt length</b>                     | 130 cm / 51.2"       | 130 cm / 51.2"       |
| <b>Colour</b>                          | black                | black                |



| 4-point pelvic harness with front-pul | Swiftly 2              |
|---------------------------------------|------------------------|
| <b>Item code</b>                      | <b>9800084711000</b>   |
| <b>Upholstered area (w x h)</b>       | 17 x 8 cm / 6.7 x 3.1" |
| <b>Rear belt length</b>               | 120 cm / 47.2"         |
| <b>Front belt length</b>              | 50 cm / 19.7"          |
| <b>Colour</b>                         | black                  |



| Crotch strap       | Swiftly   | Swiftly 2                 |
|--------------------|---|---------------------------|
| <b>Item code</b>   | <b>9800063411000</b>  | <b>9800063411000</b>      |
| <b>Belt length</b> | ca. 18 - 33 cm / 7.1 - 13"                                    | ca.18 - 33 cm / 7.1 - 13" |
| <b>Colour</b>      | black   | black                     |
| <b>Note</b>        | can only be used in combination with pelvic harness or H-belt |                           |



| Soft pelvic harness                | Swiftly              | Swiftly 2            |
|------------------------------------|----------------------|----------------------|
| <b>Item code</b>                   | <b>9800066146000</b> | <b>9800067146000</b> |
| <b>Width upholstered area rear</b> | 21 cm / 8.3"         | 26,5 cm / 10.4"      |
| <b>Belt length</b>                 | 45 cm / 17.7"        | 45 cm / 17.7"        |
| <b>Colour</b>                      | black                | black                |



| Butterfly/chest harness |  | size S               | size M               | size L               |
|-------------------------|--|----------------------|----------------------|----------------------|
| <b>Item code</b>        |  | <b>9800066166000</b> | <b>9800066176000</b> | <b>9800067176000</b> |
| <b>A</b>                |  | 18,5 cm / 7.2"       | 25 cm / 9.8"         | 25,5 cm / 10"        |
| <b>B</b>                |  | 30 cm / 11.8"        | 33,5 cm / 13.1"      | 45 cm / 17.7"        |
| <b>C</b>                |  | 16 cm / 6.3"         | 18,5 cm / 7.2"       | 23 cm / 9.1"         |
| <b>D</b>                |  | 14 cm / 5.5"         | 17 cm / 6.7"         | 22 cm / 8.7"         |
| <b>E</b>                |  | 11,5 cm / 4.5"       | 10,5 cm / 4.1"       | 11,5 cm / 4.5"       |
| <b>F</b>                |  | 18 cm / 7.1"         | 21,5 cm / 8.5"       | 25 cm / 9.8"         |
| <b>Colour</b>           |  | black                | black                | black                |



| Butterfly/chest harness, flexible |  | size S               | size M               |
|-----------------------------------|--|----------------------|----------------------|
| <b>Item code</b>                  |  | <b>9800063853000</b> | <b>9800065853000</b> |
| <b>A</b>                          |  | 12 cm / 4.7"         | 18 cm / 7.1"         |
| <b>B</b>                          |  | 21 cm / 8.3"         | 33 cm / 13"          |
| <b>C</b>                          |  | 11,5 cm / 4.5"       | 15 cm / 5.9"         |
| <b>D</b>                          |  | 9,5 cm / 3.7"        | 18 cm / 7.1"         |
| <b>E</b>                          |  | 6 cm / 3.4"          | 10 cm / 3.9"         |
| <b>F</b>                          |  | 21 cm / 8.3"         | 27 cm / 10.6"        |
| <b>Colour</b>                     |  | black                | black                |



| Seat waistcoat                  |  | size S                   | size M                   | size L                     |
|---------------------------------|--|--------------------------|--------------------------|----------------------------|
| <b>Item code</b>                |  | <b>9800066315100</b>     | <b>9800066325100</b>     | <b>9800067315100</b>       |
| <b>Upholstered area (w x l)</b> |  | 33 x 30 cm<br>13 x 11.8" | 38 x 35 cm<br>15 x 13.7" | 44 x 41 cm<br>17.3 x 16.1" |
| <b>Belt length, top</b>         |  | 65 cm / 25.6"            | 65 cm / 25.6"            | 65 cm / 25.6"              |
| <b>Colour</b>                   |  | black                    | black                    | black                      |



| Head pillow             |  | Swiftly                               | Swiftly 2                           |
|-------------------------|--|---------------------------------------|-------------------------------------|
| <b>Item code</b>        |  | <b>9800068240000</b>                  | <b>9800068370000</b>                |
| <b>Size (w x h x d)</b> |  | 26 x 13,5 x 7,5 cm<br>10.2 x 5.3 x 3" | 30,5 x 16 x 7,5 cm<br>12 x 6.3 x 3" |
| <b>Colour</b>           |  | black                                 | black                               |
| <b>Item code</b>        |  | –                                     | <b>9800068371000</b>                |
| <b>Size (w x h x d)</b> |  | –                                     | s. option black                     |
| <b>Colour</b>           |  | –                                     | red                                 |

# ACCESSORIES



| Sun and rain cover with leg blanket    | Swiftly                                      | Swiftly 2                       |
|--|--|---------------------------------|
| <b>Item code</b>                       | <b>9800066376000</b>                         | <b>9800066376000</b>            |
| <b>Shield</b> (w x h, inner dimension) | 43 x 51 cm /<br>16.9 x 20.1"                 | 43 x 51 cm /<br>16.9 x 20.1"    |
| <b>Shield bag</b> (w x h x d)          | 24 x 28 x 5 cm<br>9.4 x 11 x 2"              | 24 x 28 x 5 cm<br>9.4 x 11 x 2" |
| <b>Leg blanket</b> (w x l)             | 43 x 100 cm /<br>16.9 x 39.4"                | 43 x 100 cm /<br>16.9 x 39.4"   |
| <b>Colour</b>                          | black / light grey                           | black / light grey              |
| <b>Note</b>                            | suitable up to max. 71 cm / 28" crown height |                                 |



| Sun and rain cover                     | Swiftly                                      | Swiftly 2                       |
|--|--|---------------------------------|
| <b>Item code</b>                       | <b>9800066376100</b>                         | <b>9800066376100</b>            |
| <b>Shield</b> (w x h, inner dimension) | 43 x 51 cm /<br>16.9 x 20.1"                 | 43 x 51 cm /<br>16.9 x 20.1"    |
| <b>Shield bag</b> (w x h x d)          | 24 x 28 x 5 cm<br>9.4 x 11 x 2"              | 24 x 28 x 5 cm<br>9.4 x 11 x 2" |
| <b>Colour</b>                          | black / light grey                           | black / light grey              |
| <b>Note</b>                            | suitable up to max. 71 cm / 28" crown height |                                 |



| Rain shield, transparent for sun and rain cover | Swiftly  | Swiftly 2                     |
|---|--|-------------------------------|
| <b>Item code</b>                                | <b>9800067390000</b>   | <b>9800067390000</b>          |
| <b>Shield</b> (w x h, inner dimension)          | 43 x 51 cm /<br>16.9 x 20.1"   | 43 x 51 cm /<br>16.9 x 20.1"  |
| <b>Opening</b> (w x h)                          | 33 x 34 cm /<br>13 x 13.3"   | 33 x 34 cm /<br>13 x 13.3"    |
| <b>Leg blanket</b> (w x l)                      | 60 x 120 cm /<br>23.6 x 47.2"  | 60 x 120 cm /<br>23.6 x 47.2" |
| <b>Colour</b>                                   | transparent  | transparent                   |
| <b>Note</b>                                     | can only be used in combination with „sun and rain cover with leg blanket“ or „sun and rain cover“ |                               |





| Rain cover, transparent, with leg blanket | Swiftly | Swiftly 2 |
|---|---------|-----------|
|---|---------|-----------|

|  |  |                               |
|--|--|-------------------------------|
| <b>Item code</b>                       | <b>9800066381000</b>                         | <b>9800066381000</b>          |
| <b>Shield</b> (w x h, inner dimension) | 43 x 51 cm /<br>16.9 x 20.1"                 | 43 x 51 cm /<br>16.9 x 20.1"  |
| <b>Leg blanket</b> (w x l)             | 43 x 100 cm /<br>16.9 x 39.4"                | 43 x 100 cm /<br>16.9 x 39.4" |
| <b>Colour</b>                          | transparent                                  | transparent                   |
| <b>Note</b>                            | suitable up to max. 71 cm / 28" crown height |                               |



| Umbrella | Swiftly | Swiftly 2 |
|----------|---------|-----------|
|----------|---------|-----------|

|                   |                      |                      |
|-------------------|----------------------|----------------------|
| <b>Item code</b>  | <b>9800065940000</b> | <b>9800065940000</b> |
| <b>Diameter Ø</b> | 65 cm / 25.6"        | 65 cm / 25.6"        |
| <b>Height</b>     | 75 cm / 29.5"        | 75 cm / 29.5"        |
| <b>Colour</b>     | grey                 | grey                 |



| Sleeping bag, padded | Swiftly | Swiftly 2 |
|----------------------|---------|-----------|
|----------------------|---------|-----------|

|  |  |   |
|--|--|---|
| <b>Item code</b>                       | <b>9800066195000</b>                       | <b>9800067195000</b>                        |
| <b>Bottom part leg blanket</b> (w x l) | 35 x (68 - 80) cm<br>13.8 x (26.8 x 31.5)" | 40 x (90 - 104) cm<br>15.7 x (35.4 x 40.9)" |
| <b>Leg blanket</b> (lower leg length)  | 31 cm / 12.2"                              | 43 cm / 16.9"                               |
| <b>Colour</b>                          | black / light grey                         | black / light grey                          |


| Sleeping bag, woven fur | Swiftly | Swiftly 2 |
|-------------------------|---------|-----------|
|-------------------------|---------|-----------|

|  |  |   |
|--|--|---|
| <b>Item code</b>                       | <b>9800066185000</b>                       | <b>9800067185000</b>                        |
| <b>Bottom part leg blanket</b> (w x l) | 35 x (68 - 80) cm<br>13.8 x (26.8 x 31.5)" | 40 x (90 - 104) cm<br>15.7 x (35.4 x 40.9)" |
| <b>Leg blanket</b> (lower leg length)  | 31 cm / 12.2"                              | 43 cm / 16.9"                               |
| <b>Colour</b>                          | black / light grey                         | black / light grey                          |

# ACCESSORIES



| Basket                  | Swiftly                                 |
|-------------------------|---|
| <b>Item code</b>        | <b>9800068130000</b>                    |
| <b>Size (w x h x l)</b> | 26 x 17 x 30 cm /<br>10.2 x 6.7 x 11.8" |
| <b>Volume</b>           | 13 l                                    |
| <b>Loading capacity</b> | 3 kg / 6.6 lbs                          |
| <b>Colour</b>           | black                                   |



| Basket                  | Swiftly 2                             |
|-------------------------|---------------------------------------|
| <b>Item code</b>        | <b>9800068350000</b>                  |
| <b>Size (w x h x l)</b> | 33 x 18,5 x 30 cm<br>13 x 7.3 x 11.8" |
| <b>Volume</b>           | 15 l                                  |
| <b>Loading capacity</b> | 3 kg / lbs                            |
| <b>Colour</b>           | black                                 |



| Pneumatic wheels, upgrade kit           | Swiftly              |
|---|----------------------|
| <b>Item code</b>                        | <b>9800066470013</b> |
| <b>Wheel size swivel wheels (front)</b> | 7,5"                 |
| <b>Wheel size (rear)</b>                | 10"                  |



| Pneumatic wheels, upgrade kit           | Swiftly 2            |
|---|----------------------|
| <b>Item code</b>                        | <b>9800066710000</b> |
| <b>Wheel size swivel wheels (front)</b> | 7,5"                 |
| <b>Wheel size (rear)</b>                | 10"                  |





| Tip assist                | Swiftly               | Swiftly 2              |
|---------------------------|-----------------------|------------------------|
| <b>Item code</b>          | <b>9800068291000</b>  | <b>9800068291000</b>   |
| <b>Tread area (w x l)</b> | 4 x 8 cm / 1.6 x 3.1" | 4 x 8 cm / 1.6 x 3.1 " |
| <b>Colour</b>             | black                 | black                  |



Member of the International  
Support Association for the  
Rehabilitation of Children and Youths e.V.

---



Thomas Hilfen für Körperbehinderte  
GmbH & Co. Medico KG

Walkmühlenstraße 1 · 27432 Bremervörde · Germany  
Phone: +49 (0) 4761 8860 · Fax: +49 (0) 4761 886-19  
info@thomashilfen.de · www.thomashilfen.com

**USA and Canada:**

Thomashilfen North America  
14221 29th St E, #101 · Sumner, WA 98390 · USA  
phone: 866 870 2122 (toll free in U.S.)  
info@thomashilfen.com · www.thomashilfen.us