

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

**Step 1 – Seat Sizes (Required options are indicated below by † symbol)**

	ThevoTwist seat size 1			ThevoTwist seat size 2			HCPCS	Description
	£ Blue	£ Red	\$420	£ Blue	£ Red	\$420		
†Seat	668821	668811	\$420	678821	678811	\$420	E2292	Individual body adapted, integrated seat/back with MIS-Microstimulation®; 145° recline gas spring operated; growth adjustable seat upholstery; articulated footrest height and angle adjustable; height, depth angle adjustable swing away head rest bracket; front/rear facing seat; ventilated/washable upholstery in 2 colors
†Back	668822	668812	\$518	678822	678812	\$518	E2293	
†Recline	668823	668813	\$415	678823	678813	\$415	E1014	
†Hip sup.	668824	668814	\$230	678824	678814	\$230	E0956	
†Calf panel	6635	6635	\$45	6735	6735	\$45	E0995	
†Foot plate, flip up	668825	668815	\$90	678825	678815	\$90	K0040	
†Legrest, elevating	668826	668816	\$200	678826	678816	\$200	E0990	



6387/ 6587 £†Standard      6598 £†Neck support      6599 £†Shell

\$285	E0955	†Headrest is a required option!
\$170	E1028	†Headrest bracket is a required option!

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Seat to Head	Max Height	Hip Angle
Size 1	7 – 11½"	6¾ – 12"	7¾ – 11"	18 ¼ – 26 ¾"	49"	90 - 133°
Size 2	8½ – 14½"	9¼ – 15"	10 – 13½"	25 ¼ – 30 ¾"	60"	90 - 145°

**Step 2 – Frame Options (Required options are indicated below by † symbol)**

**A-chassis** - Folding aluminum frame, tilt-in-space, dual spring suspension, cargo basket, tilt-in-space lever at push bar

**Q-chassis** - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas spring, lockable rear wheels



	HCPCS	SIZE 1	SIZE 2
A-chassis	E1234	£ 6202    \$2295	£ 6402    \$2295
†Castor pin lock, pair	K0073	67830    \$96	67830    \$96
†Bus tie-down loop kit (*)	K0108	67839    \$205	67839    \$205
Pneumatic wheels, full set	K0070	£ 6646    \$249	£ 6646    \$249
Frame padding, black	K0108	£ 6622    \$115	£ 6622    \$115
Frame padding, red/ black	K0108	£ 66221    \$115	£ 66221    \$115
Tip Assist	K0108	£ 6750    \$89	£ 6750    \$89
Under seat storage basket	K0108	NC	NC

	HCPCS	SIZE 1 & 2
Q-chassis, hi/low indoor base, pedal op.	E1399	£ 6666    \$1999
Push bar for Q-chassis	K0108	£ 6660    \$298

**Q-chassis Technical Specs \*\***

Capacity	88 lbs
Seat tilt	-7 to +31°
Floor to seat	11 - 28"
Seat direction	front facing

**ThevoTwist Technical Specs \*\***

	SIZE 1	SIZE 2
Seat capacity	77 lbs	88 lbs
Frame capacity	99 lbs	110 lbs
Unit weight	36 lbs	42 lbs
Max seat tilt	31°	
Seat suspension, lockable	Yes	
Upholstery	Trevira-CS	
Foot angle	78° - 100°	
Size of foot rest	11½" x 7⅞"	13" x 8¼"
Floor to handle	30 – 46"	32 – 47"
Floor to seat	19"	19½"
Max width	23"	27"
Fold length	32"	35½"
Fold height	18½"	18½"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7½" solid	
Rear wheels	10" solid	

**Z-chassis** - indoor base/feeder/therapy frame, manual height/ tilt adjustable, lockable rear wheels



	HCPCS	SIZE 1 & 2
Z-chassis	E1399	£ 6691    \$899

**Z-chassis Technical Specs \*\***

Capacity	88 lbs
Seat tilt	-15 to +35°
Floor to seat	15 - 23"
Seat direction	front facing

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

\*\*More details in Technical Manual

**Step 3 – Positioning & General Options (Required options are indicated below by † symbol)**

<b>Lower Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle		E0952	Pair		£ 6648		\$35
Ankle Huggers, Bodypoint®		K0108	Pair		£ 213151		\$205
Footrest pad, height 1¾"		K0108	Each	£ 10622		£ 10023	\$99

<b>Mid Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Abduction block, padded		E0957	Each	£ 6620		£ 6720	\$160
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420		£ 8470	\$95
	4 point, padded, SR buckle	E0978	Each	£ 8421		£ 8471	\$135
	Soft pelvic harness	E0978	Each	£ 6614		£ 6714	\$175
Crotch strap (fits 8420,8470,8421,8471)		E0978	Each		£ 6341		\$33
Tray, height & depth adjustable, clear w/picture window		E0950	Each	£ 66281		£ 67281	\$190
Mounting hardware tray – †required when ordering 6628/6728		E1028	Each	£ 66282		£ 67282	\$225
Arm rest, height, depth and angle adjustable		E0973	Pair	£ 8431		£ 8478	\$375

<b>Upper Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385		£ 6585	\$189
	Vest style	E0960	Each	£ 6631	£ 6632	£ 6731	\$189
	H-style, posture support, small, Bodypoint®	E0960	Each		£ 214806		\$199
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		£ 214270		\$199
	PivotFit shoulder harness, Bodypoint®	E0960	Each		£ 214250		\$189
Trunk supports	Lateral, Dynamic, MIS-Microstimulation®			£ 6383		£ 6583	\$255
	Rigid lateral padded, left & right, height adjustable	E0956	Pair	£ 6377		£ 6577	\$245
	Flexible lateral, padded, w/chest strap, height adj.	E0956	Pair		£ 6596		\$295

<b>General Accessories</b>		HCPCS	Unit	Small	Medium	Large	Price
Grab bar, padded, black		K0108	Each	£ 6621		£ 6721	\$125
Vent tray		E1029	Each		£ 110004		\$450
Oxygen tank holder, E-tank		E2208	Each		£ 665012		\$210
Straps	Secure folded seat/frame	K0108	Each		£ 667011		\$30
	Wrist "runaway" tether	K0108	Each		£ 667111		\$30
Weather protection	Canopy, sun/rain with window, black	K0108	Each		£ 663761		\$325
	Lower body windbreaker, black for 663761	K0108	Each		£ 661852		\$99
	Canopy, sun/rain with window, grey	K0108	Each		£ 637811		\$325
	Lower body windbreaker grey for 637811	K0108	Each		£ 637812		\$99
	Canopy, rain only, clear	K0108	Each		£ 6638		\$245
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		£ 6739		\$125
	Rain cape, grey	K0108	Each	£ 6340		£ 6540	\$163
	Sleeping bag, summer weight, black	K0108	Each	£ 6619		£ 6719	\$299
	Sleeping bag, winter weight, black	K0108	Each	£ 6618		£ 6718	\$299
	Umbrella, angle adjustable	K0108	Each		£ 6594		\$95

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

\*\*More details in Technical Manual

