

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	

Step 1 – Seat Sizes (Required options are indicated below by † symbol)



	£ tRide seat size 1		£ tRide seat size 2		HCPCS	Description
†Seat	622301	\$360	642301	\$360	E2294	Solid contoured, integrated seat/back, 140° recline, growth adjustable without tools; articulated footrest height and angle adjustable; front facing seat; height and angle adjustable padded arm supports; height, depth and angle adjustable swing-away bracket with head rest neck support, tilt-in-space lever at seat unit, 0 degrees tilt stop, ventilated/ washable upholstery in 3 colors
†Back	622302	\$350	642302	\$350	E2293	
†Recline	622303	\$200	642303	\$200	E1014	
†Hip sup.	622304	\$200	642304	\$200	E0956	
†Bracket	622305	\$125	642305	\$125	E1028	
†Head rest	622306	\$245	642306	\$245	E0955	
†Calf panel	6635	\$45	6735	\$45	E0995	

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Back Height#	Max Height	Hip Angle
Size 1	6¾ – 11¾"	6¾ – 11"	7¾ – 11¾"	10¼ – 16½"	50"	90° - 140°
Size 2	9¾ – 14¼"	9½ – 13¾"	10¼ – 13¾"	12½ – 20½"	60"	90° - 140°

Back Height without headrest

†Color Options



Grey/ pink
Size 1 £ 6290
Size 2 £ 6490



Black/ grey
Size 1 £ 6291
Size 2 £ 6491



Grey/ green
Size 1 £ 6292
Size 2 £ 6492

Step 2 – Frame Options (Required options are indicated below by † symbol)

t-chassis - Folding aluminum frame, 35° posterior tilt-in-space, 0° anterior tilt stop, cargo basket (Required options are indicated below by † symbol)



Q-chassis - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas spring, lockable rear wheels (Required options are indicated below by † symbol)



	HCPCS	SIZE 1		SIZE 2	
t-chassis	E1234	£ 622307	\$2025	£ 642307	\$2025
†Castor pin lock, pair	K0073	67830	\$96	67830	\$96
†Bus tie-down loop kit (*)	K0108	67839	\$205	67839	\$205
Pneumatic wheels, full set	K0070	£ 6646	\$249	£ 6646	\$249
Frame padding, black	K0108	£ 6622	\$115	£ 6622	\$115
Frame padding, red/ black	K0108	£ 66221	\$115	£ 66221	\$115
Tip Assist/ Curb tipper	K0108	£ 6750	\$89	£ 6750	\$89
Under seat storage basket	K0108		NC		NC

	HCPCS	SIZE 1		SIZE 2	
Q-chassis, indoor hi/low	E1399	£ 6666		\$1999	
†tRide adaptor	K0108	£ 6295	\$130	£ 6495	\$130
Push bar for Q-chassis	K0108	£ 6660		\$288	

tRide Technical Specs **	SIZE 1	SIZE 2
Seat capacity	77 lbs	88 lbs
Frame capacity	99 lbs	110 lbs
Unit weight	36 lbs	42 lbs
Max seat tilt	35°	
Upholstery	Trevira-CS	
Height adjustment arm support	3¾ – 5½"	4¾ – 7½"
Foot angle	78° - 100°	
Size of foot rest	11½" x 7⅞"	13" x 8½"
Floor to handle	30 – 46⅞"	33 – 47½"
Floor to seat	18⅞"	18½"
Max width	23½"	26½"
Fold length	32"	35½"
Fold height	18½"	18½"
Brakes	foot operated	
Seat direction	front facing	
Front wheels	7½" solid	
Rear wheels	10" solid	

Q-chassis Technical Specs **	
Capacity	88 lbs
Seat tilt	-7 to +31°
Floor to seat	11 - 28"
Seat direction	front facing



* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement. ExoMotion LLC dba Thomashilfen www.thomashilfen.us Effective Date January 1st 2022

Step 3 – Positioning & General Options (Required options are indicated below by † symbol)

Lower Body Positioning		HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle		E0952	Pair		£ 6648		\$35
Ankle Huggers, Bodypoint®		K0108	Pair		£ 213151		\$205
Footrest pad, height 1¾"		K0108	Each	£ 10622		£ 10023	\$99

Mid Body Positioning

Abduction block, padded		E0957	Each	£ 6620		£ 6720	\$160
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420		£ 8470	\$95
	4 point, padded, SR buckle	E0978	Each	£ 8421		£ 8471	\$135
	Soft pelvic harness / leg harness	E0978	Each	£ 6614		£ 6714	\$175
Crotch strap (fits 6629,8420,8470,8421,8471)		E0978	Each		£ 6341		\$33
Tray, height adjustable, clear with picture window		E0950	Each	£ 66281		£ 67281	\$190
Mounting hardware tray – †required when ordering 6628/6728		E1028	Each	£ 66282		£ 67282	\$225

Upper Body Positioning

Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385		£ 6585	\$189
	Vest style	E0960	Each	£ 6631	£ 6632	£ 6731	\$189
	H-style, posture support, small, Bodypoint®	E0960	Each		£ 214806		\$199
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		£ 214270		\$199
	PivotFit shoulder harness, Bodypoint®	E0960	Each		£ 214250		\$189
Trunk supports	Rigid lateral, left & right, height adjustable	E0956	Pair	£ 6634		£ 6734	\$245
	Flexible lateral, l&r, w/ chest strap, height adj.	E0956	Pair		£ 6633		\$295
Head support	Fitting, ball joint (to fit after market head rest)	K0108	Each		£ 6268		\$90

General Accessories

Grab bar, padded, black		K0108	Each	£ 6621		£ 6721	\$125
Vent tray		E1029	Each		£ 110004		\$450
Oxygen tank holder, E-tank		E2208	Each		£ 665012		\$210
Straps	Secure folded seat/frame	K0108	Each		£ 667011		\$30
	Wrist "runaway" tether	K0108	Each		£ 667111		\$30
Weather protection	Canopy, sun/rain with window, black	K0108	Each		£ 663761		\$325
	Lower body windbreaker, black for 663761	K0108	Each		£ 661852		\$99
	Canopy, rain only, clear	K0108	Each		£ 6638		\$245
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		£ 6739		\$125
	Rain cape, grey	K0108	Each	£ 6340		£ 6540	\$163
	Sleeping bag, summer weight, black	K0108	Each	£ 6619		£ 6719	\$299
	Sleeping bag, winter weight, black	K0108	Each	£ 6618		£ 6718	\$299
Umbrella, angle adjustable		K0108	Each		£ 6594		\$95

*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

**More details in Technical Manual

